




Dolphin

Large Print Keyboard

Complete
list of 'quick
buttons' included

The quickest
way to get
started with
SuperNova &
GuideConnect



© Wacker Chemie AG

Designed for people with sight loss

The Dolphin Keyboard features large print, high contrast keys that are easy to see. It is available in both the UK and US International editions.

Main Benefits:

- High contrast, large print keys
- 18 easy to locate, tactile 'Quick Buttons'
- Fewer hot keys to remember
- Ideal for both novice and experienced keyboard users
- Get started with SuperNova & GuideConnect more quickly

For more information:

01905 754 577



Dolphin
Making a difference



Dolphin

Large Print Keyboard

The Dolphin Keyboard features large print, high contrast keys that are easy to see and also includes 18 Quick Buttons for speedy access to your Dolphin software's most useful features.

These buttons are large and of varied shapes, making them tactile and easy to locate by touch.




What you need

The Dolphin Keyboard is recommended for use with SuperNova version 13.50 or later, or GuideConnect version 1.01 or later.



















Simply plug in your Dolphin Keyboard using the USB cable attached and get started straight away!



















Image of the letter 'A' on standard keyboard compared directly with the same letter on the Dolphin Keyboard:



Warranty

Your Dolphin Keyboard comes with a 12 month warranty against manufacturing defects, which are very rare. If you suspect a problem, please contact us immediately stating where and when you purchased it. We will happily replace any keyboard with a genuine manufacturing fault. We will not replace keyboards which have been misused or damaged.

	Dolphin logo	Launch SuperNova's control panel
	Question mark	Open help
	Podcast symbol	Launch pad (podcasts, radio, etc.)
	Squares	Colour changer on/off
	Stripes	Focus highlighting on/off
	Magnifying glass	Magnification on/off
	Plus sign	Increase magnification (zoom in)
	Minus sign	Decrease magnification (zoom out)
	Monitor	Cycle through magnification modes
	Scanner	Open scan and read menu
	Triangle (play)	Read from here (press escape to stop)
	Mouth	Speech on/off
	Tortoise	Speech slower
	Hare	Speech faster
	Envelope	Open email application
	Globe	Open web browser
	Right arrow	Move forward (web browsing)
	Left arrow	Move back (web browsing)

	Dolphin logo	Launch GuideConnect or Open Main Menu if GuideConnect already open
	Question mark	Open help
	Podcast symbol	Open entertainment menu
	Squares	Cycle colour scheme
	Stripes	Open actions menu
	Magnifying glass	Not supported
	Plus sign	Increase magnification (zoom in)
	Minus sign	Decrease magnification (zoom out)
	Monitor	Not supported
	Scanner	Open scanner menu
	Triangle (play)	Continuous read
	Mouth	Voice on/off
	Tortoise	Decrease voice speed
	Hare	Increase voice speed
	Envelope	Open email menu
	Globe	Opens website menu
	Right arrow	Forward a page in web browser
	Left arrow	Back a page in web browser